

## Sampling Dipper with integral handle



**Model : YR1064**

**Sampling Dipper with integral handle**

**Manufacturer : Kalstein**

**Price : \$ 0.00**

**Description :**

For any laboratory, it is fundamental that the tests they make will not be affected by any outer element.

The sample dipping process is an important step in any laboratory. Which means that it is necessary to be careful when it comes to handling the samples.

### Sampling dipper

A sampling dipper is a very important tool inside a laboratory used to collect and dispense any kind of liquids in a safely way. They are also useful for transferring hot or cold reactive substances. These dippers are designed for sterile sample collection and dispatch so there is no risk of cross contamination.

### Water sample dippers

These lightweight dippers are easy to use and will not corrode or contaminate solutions. They have polypropylene caps, which is a thermoplastic polymer resistant to many chemical solvents, bases and acids. Polypropylene's melting point is 130 degrees Celsius. These dippers have a polystyrene container, which is a suitable material for all kinds of thermoformed packaging, providing the necessary safety for any containers.

### Stainless steel sampling dipper

As you may know by now, these are not stainless steel models. One of the great features of these sampling dippers is that they come with an integrated handle, which makes it very practical for the everyday use in a laboratory.

You should know that they are individually wrapped and processed with an Ethylene Oxide Sterilization. In addition, they have a "flow seal" screw cap, which are leak-proof tested in accordance with EN14245 Annex D and BS5213.

The sampling dipper's containers can come in a clear or a blue presentation. The blue containers are highly visible if accidentally dropped into a liquid sample. There are 30ml models of clear and blue containers with a polypropylene cap. You should know that there is a 100ml model of

## Sampling Dipper with integral handle

clear and blue containers with metal cap. Finally, a 250ml model of clear and blue containers with metal cap. These sampling dippers come distributed as 50 units per case.

Designed for sterile sample collection and dispatch to avoid the risk of cross contamination.

Individually wrapped and EO sterility.

Blue containers are highly visible if accidentally dropped into a liquid sample.

"Flow seal" screw cap, leakproof-tested in accordance with EN14245 Annex D and BS5213.

Parameters by Models

## Sampling Dipper with integral handle

Sampling Dipper with integral handle				
Model	Description	Sterility	Inner packing	Case Qty.
YR1064	Sampling Dipper with integral handle, 30mL, PP cap, PS, Clear	EO	1/flow pack	50
YR1065	Sampling Dipper with integral handle, 100mL, Metal cap, PS, Clear	EO	1/flow pack	50
YR1066	Sampling Dipper with integral handle, 250mL, Metal cap, PS, Clear	EO	1/flow pack	50

## Sampling Dipper with integral handle

### Options available for Sampling Dipper with integral handle :

#### Variation

- (YR1064) Capacity: 30ml / Cap: PP,
- (YR1065) Capacity: 100ml / Cap: Metal (+\$ 23.37),
- (YR1066) Capacity: 250ml / Cap: Metal (+\$ 43.81).