

## Specimen container, Yellow screw cap assembled, 60ml



**Model : YR1069**

**Specimen container, Yellow screw cap assembled, 60ml**

**Manufacturer : Kalstein**

**Price : \$ 108.09**

### **Description :**

For any lab, it is fundamental that the tests they make is safe. The handled specimen has to go through a very delicate process, before and after the tests. It is necessary to be extremely careful, and one way to accomplish that is by getting the right specimen containers.

The transportation of specimen and delicate substances needs to be in a proper way to avoid any leaks that can provoke the loss of important matter, so it is fundamental to have the proper container in order to receive good turnouts. The 60ml specimen container with yellow screw cap is the right choice for collection and transportation of delicate substances.

### Specimen containers

Polypropylene is very common in industrial applications due to its unique properties and ability to adapt, especially in the packaging area. This is a virgin polypropylene container, which makes it the ideal decision for the collection and transportation of any kind of specimen. One of the most important features of this container is that it does not break. Besides, polypropylene is a "thermoplastic", which means that this container can reach a temperature of 130 degrees Celsius (being this the melting point).

### Non-cytotoxic

Due to the goals of this product, this container is not toxic to cells, making it the right choice for collection and transportation of specimen. In other matter, the dual threaded caps reinforces the seal, so that any type of leaks can be avoided, giving the container a great resistance seal, in that way it can be preserved and moved around in harsh conditions.

### Specimen containers with lids

In any laboratory, the results of the tests rely on the way the containers are organized and stored. Therefore, coding and storage is very important. These ones have a variety of colored caps, which make them easier to classify and put in storage, being very appropriate for coding in laboratories. These specimen containers with red security seal are easy stackable and compact so that storage will not be a problem.

There are four types of these specimen containers, two being with no label, one of them sterile and the other non-sterile. Moreover, there is the other

## **Specimen container, Yellow screw cap assembled, 60ml**

couple of containers with label, one of them sterile and the other non-sterile

Manufactured from virgin polypropylene. All materials are non-cytotoxic.

Ideal for specimen collection and transportation in harsh conditions.

Dual threaded caps for a leak resistant seal.

Available in different cap colors for visibility or coding.

Containers are easy stackable and compact for easy storage.

Parameters by Models

**Specimen container, Yellow screw cap assembled, 60ml****Specimen container, Yellow screw cap assembled**

Model	Description	Sterility	Inner packing	Case Qty.
YR1069	Specimen container, Yellow screw cap, Assembled, PP, No label	N.S.	Bulk	500
YR1070	Specimen container, Yellow screw cap, Assembled, PP, With label	EO	25/zip bag	500

## Specimen container, Yellow screw cap assembled, 60ml

**Options available for Specimen container, Yellow screw cap assembled, 60ml :**

Variation

(YR1069) No label,

(YR1070) With label (+\$ 32.13).