

Mould Incubator



Model : YR0554

Mould Incubator

Manufacturer : Kalstein

Price : \$ 3,200.00

Description :

Incubators are an essential dispositive in a laboratory when it comes to cellular biology experiments, microbiology, and cellular cultures. The perfect incubator maintains the right temperature, humidity and other factors involved in experiments that imitate the living process. Kalstein products offers you the ideal mould incubator so that you can rely on the storage of any kind of culture you want to preserve.

These mould incubator have terrific features that guarantees the safety and quality of the cultures you want to maintain. They very common in the fields of pharmacy, chemical industry, environmental protection, public health, agriculture and many more; they are special equipment to cultivate organisms for scientific research and examining.

Get to know some of this device's great features. First, the body and chamber adopt the round angle structure; this makes it even easier to clean. This mould incubator also has a power socket in the chamber and a UV light for the sterilization process. About the safety of your samples, the device has a double door design and a magnetic door sealing. In addition, it has a cooling / heating automatic control. In addition, this mould incubator runs with the International famous compressor, CFC free refrigerant high efficient and low energy consumption.

Features:

- The body and the chamber adopt the round angle structure, easy to clean.
- There is power socket in the chamber, and UV lamp for sterilization.
- Double door design and magnetic door sealing.
- Cooling /heating automatic control.
- International famous compressor, CFC free refrigerant, high efficient and low energy consumption.

Optional accessories:

- Printer (For: YR0559, YR0560, YR0561, YR0562, YR0563).
- Independent temperature protector.

Standard accessory:

Mould Incubator

- Humidifying plate (For: YR0554, YR0555, YR0556, YR0557, YR0558).
- Humidifier, micro printer (For: YR0564, YR0565, YR0566, YR0567, YR0568).

Parameters by Models

Mould Incubator

Mould Incubator

Model:	YR0554	YR0555	YR0556	YR0557	YR0558
Capacity:	150L	200L	250L	300L	400L
Temp.range:	5~50°C	5~50°C	5~50°C	5~50°C	5~50°C
Temp.Fluctuation:	±1 °C	±1 °C	±1 °C	±1 °C	±1 °C
Shelves:	3 pcs/ adjustable	3 pcs/ adjustable	3 pcs/ adjustable	3 pcs/ adjustable	3 pcs/ adjustable
Consumption:	250W	300W	350W	400W	500W
Power supply:	110/220V±10%, 50/60Hz 110/220V±10%, 50/60Hz 110/220V±10%, 50/60Hz 110/220V±10%, 50/60Hz 110/220V±10%, 50/60Hz				
Standard accessory:	Humidifying plate				
Internal size(mm):	450*420*850	450*420*850	480*490*1070	520*500*1170	580*540*1270
External size(mm):	590*620*1380	590*640*1520	620*680*1590	660*690*1690	720*740*1790
Packing size(mm):	750*780*1560	750*800*1710	780*840*1780	820*850*1880	880*900*1980
Gross weight(kg):	120	138	147	157	170

Mould Incubator

Mould Incubator					
Model:	YR0559, YR0564	YR0560, YR0565	YR0561, YR0566	YR0562, YR0567	YR0563, YR0568
Capacity:	150L	200L	250L	300L	400L
Temp. Range:	0~65°C	0~65°C	0~65°C	0~65°C	0~65°C
Temp. Fluctuation:	±1 °C	±1 °C	±1 °C	±1 °C	±1 °C
Shelves:	2 pcs/adjustable	2 pcs/adjustable	2 pcs/adjustable	2 pcs/adjustable	2 pcs/adjustable
Consumption:	300W	350W	430W	480W	750W
Power supply:	110/220V±10%, 50/60Hz 110/220V±10%, 50/60Hz 110/220V±10%, 50/60Hz 110/220V±10%, 50/60Hz 110/220V±10%, 50/60Hz				
Standard accessory:	Humidifier, micro printer (For: YR0564, YR0565, YR0566, YR0567, YR0568)				
Internal size (mm):	450*420*850	450*450*1000	480*490*1070	520*500*1170	580*540*1270
External size (mm):	630*630*1410	630*630*1570	680*630*1640	700*670*1750	730*730*1830
Packing size (mm):	770*770*1590	770*770*1750	820*770*1810	840*800*1920	870*870*2000
Gross weight (kg):	120	138	147	157	170

Mould Incubator

Options available for Mould Incubator :

Variation

- (YR0558) 400L / With humidifying plate (+\$ 1,300.00),
- (YR0559) 150L / With humidifier (+\$ 1,450.00),
- (YR0560) 200L / With humidifier (+\$ 1,900.00),
- (YR0561) 250L / With humidifier (+\$ 2,350.00),
- (YR0562) 300L / With humidifier (+\$ 2,850.00),
- (YR0563) 400L / With humidifier (+\$ 2,850.00),
- (YR0564) 150L / With micro printer (+\$ 2,200.00),
- (YR0554) 150L / With humidifying plate,
- (YR0565) 200L / With micro printer (+\$ 2,550.00),
- (YR0555) 200L / With humidifying plate (+\$ 150.00),
- (YR0566) 250L / With micro printer (+\$ 3,000.00),
- (YR0556) 250L / With humidifying plate (+\$ 350.00),
- (YR0567) 300L / With micro printer (+\$ 3,500.00),
- (YR0557) 300L / With humidifying plate (+\$ 750.00),
- (YR0568) 400L / With micro printer (+\$ 4,450.00).