

## Transfer pipette, 5ml Capacity-Graduated to 1ml - Large Bulb



**Model : YR1119**

**Transfer pipette, 5ml Capacity-Graduated to 1ml - Large Bulb**

**Manufacturer : Kalstein**

**Price : \$ 82.38**

**Description :**

One of the most common tools and devices inside any lab is the pipette. This is a lab tool, commonly used in chemistry, biology and medicine to transport a measured volume of liquid, often as a media dispenser. There are several designs and purposes with different levels of accuracy and precision.

### Transfer pipette 5ml

This kind is the most basic pipettes in the area. You can use it for rough measurements only; however, it is still very important that you follow a standard pipetting technique. This is a single-piece pipette and its bulb can serve as the liquid-holding chamber.

### Transfer pipettes function

About the functions, it is recommended that the user aspirates the liquid at a 90 degree angle, dispense at a 45 degree angle and touch off to make sure all the liquid is dispensed. The uses ranges from general purpose to sensitive lab research.

### Transfer pipette 1ml

The pipette we are presenting to you has a capacity of 5ml and is 1ml graduated, with a graduated increase of 0.25ml. The total length of the transfer pipette is 160mm. The stem diameter is 5.0mm. This is a large bulb transfer pipette so the bulb draw is 3.5 ml. Finally, the dropper/ml is 25.

About the packaging, the transfer pipettes can come as 500 unit per box and 3000 per case, this category is the non-sterile type. There is also the option of packaging 1 per peel pack and 500 per case, 5 per peel pack and 500 per case, 10 per peel pack, 1500 per case, and 20 per peel pack and 1500 per case. These last four come with Ethylene Oxide Sterilization process.

### Parameters by Models

## Transfer pipette, 5ml Capacity-Graduated to 1ml - Large Bulb

Transfer pipette, 5ml Capacity-Graduated to 1ml - Large Bulb

Model	Description
YR1119	500/Box,3000/Case,Nonsterile
YR1120	1/Peel Pack,500/Case,Sterile
YR1121	5/Peel Pack,500/Case,Sterile

## Transfer pipette, 5ml Capacity-Graduated to 1ml - Large Bulb

### Options available for Transfer pipette, 5ml Capacity-Graduated to 1ml - Large Bulb :

#### Variation

(YR1120) Peel Pack: 1 / Case: 500 / Sterile (-\$ 38.56),  
(YR1121) Peel Pack: 5 / Case: 500 / Sterile (-\$ 44.40),  
(YR1119) Box: 500 / Case: 3000 / Non sterile.