

Transfer pipette, 5ml Capacity-Graduated to 2ml - Blood Bank



Model : YR1125

Transfer pipette, 5ml Capacity-Graduated to 2ml - Blood Bank

Manufacturer : Kalstein

Price : \$ 91.14

Description :

When it comes to lab material, there is an infinity of items necessary in the area. Without them, testing samples and working in any medical environment would not be impossible. Kalstein products is offering you the possibility of having the high-end tools and products for your laboratory, so that you can rely and trust that the work you are doing, is the best.

One of the most common tools and devices inside any lab is the pipette. This is a lab tool, commonly used in chemistry, biology and medicine to transport a measured volume of liquid, often as a media dispenser. There are several designs and purposes with different levels of accuracy and precision.

Transfer pipette

This kind is the most basic pipettes in the area. You can use it for rough measurements only; however, it is still very important that you follow a standard pipetting technique. This is a single piece of plastic pipette, and its bulb serves as the liquid-holding chamber.

Transfer pipettes function

About the functions, it is recommended that the user aspirates the liquid at a 90 degree angle, dispense at a 45 degree angle and touch off to make sure all the liquid is dispensed. The uses ranges from general purpose to sensitive lab research.

Transfer pipette 5ml

The pipette we are presenting to you has a capacity of 5ml and is 2ml graduated with a graduated increasement of 0.5ml. The total length of the transfer pipette is 155mm. The stem diameter is 6.7mm. This is a large bulb blood bank transfer pipette so the bulb draw is 2.4 ml. Finally, the dropper/ml is 22.

About the packaging, the transfer pipettes can come as 500 unit per box and 3000 per case, this category is the non-sterile type. There is also the option of packaging 1 per peel pack and 500 per case, 5 per peel pack and 500 per case, 10 per peel pack, 1500 per case, and 20 per peel pack and 1500 per case. These last four come with an Ethylene Oxide Sterilization.

Transfer pipette, 5ml Capacity-Graduated to 2ml - Blood Bank

Parameters by Models

Transfer pipette, 5ml Capacity-Graduated to 2ml - Blood Bank**Transfer pipette, 5ml Capacity-Graduated to 2ml - Blood Bank**

Model	Description
YR1125	500/Box,3000/Case,Nonsterile
YR1126	1/Peel Pack,500/Case,Sterile
YR1127	5/Peel Pack,500/Case,Sterile

Transfer pipette, 5ml Capacity-Graduated to 2ml - Blood Bank

Options available for Transfer pipette, 5ml Capacity-Graduated to 2ml - Blood Bank :

Variation

- (YR1127) Peel Pack: 5 / Case: 500 / Sterile (-\$ 53.17),
- (YR1125) Box: 500 / Case: 3000 / Non sterile,
- (YR1126) Peel Pack: 1 / Case: 500 / Sterile (-\$ 47.32).