


Pipette Tips



Model : YR1241

Pipette Tips

Manufacturer : Kalstein

Price : \$ 315.81

Description :

Is very common to find pipettes in any chemistry, biology or medical laboratory, since they are a common laboratory tool to transfer a measured volume of liquid. Of course, they come in different designs and different purposes. Regardless of the type of pipette, you choose you need to have a proper pipette tip in order to have accurate results during the testing process. Many pipette types work by creating a partial vacuum above the liquid-holding chamber and selectively releasing this vacuum to draw up and dispense liquid.

Pipette tips functions

Since liquids are never drawn directly into the pipette's shaft, you need to have the right pipette's tips, which are disposable plastic objects attached to the shaft. Because of their importance, you need to choose the right pipette tips in order to have good and reliable results. That is why you need to select the right ones and handle them carefully.

Pipette tips 1000ul

In any laboratory, liquid handling is a very common practice when it comes to sample procedures. In order to transport an exact amount of liquid, you must have pipettes. That is why you should know that good pipetting results rely on the combination of a good quality tip and the right pipette. The pipette tips we are offering you, will provide maximum accuracy and precision.

Pipette tips yellow

These prime quality virgin polypropylene tips' design fits a wide variety of pipettors. You should also know, that polypropylene is a thermoplastic polymer resistant to many chemical solvents, bases and acids. This particular feature makes these boxes the most suitable for keeping you lab ware safe. Another important factor is that polypropylene's melting point is 130 degrees Celsius. These tips have the right shape to ensure precise measurements, being a high-end product. In addition, these tips come packaged in re-sealable bags.

We provide you with different types, colors and measurements pipette tips. We have them available from 5 – 10 ul, 1 – 200 ul, 10 – 300 ul, 5 – 200 ul to 1000 ul. Among the types you can find universal Gilson, Eppendorf with

Pipette Tips

collar, crystal graduated, beveled, universal Eppendorf, Eppendorf crystal, MLA, Oxford and extra-long. About the colors, you can find them in three different options natural, yellow or blue. These pipette tips come packed from 30000 per case to 10000 per case.

- Tips designed to fit a wide variety of pipettors. Manufactured from prime quality virgin polypropylene and molded to exact standards to ensure precise measurements.
- Packaged in resealable bags.

Parameters by Models

Pipette Tips

Pipette Tips				
Model	Description	Sterility	Inner packing	Case Qty.
YR1241	Pipette tip, 5 -10 ul,Gilson crystal, graduated, Natural	N.S.	1000	30000
YR1242	Pipette tip, 1 -200 ul,Gilson graduated, bevelled.Natural	N.S.	1000	30000
YR1243	Pipette tip, 1 -200 ul,Gilson graduated, bevelled.Yellow	N.S.	1000	30000
YR1244	Pipette tip, 5 -200ul,Universal Gilson, Yellow	N.S.	1000	30000
YR1245	Pipette tip, 100 -1000 ul,Universal Gilson,graduated, bevelled	N.S.	500	10000
YR1246	Pipette tip, 0,5 -10 ul,Eppendorf crystal,Natural	N.S.	1000	30000
YR1247	Pipette tip, 5 -200ul,Eppendorf with collar, Yellow	N.S.	1000	30000
YR1248	Pipette tip, 10 -300ul,Universal Eppendorf graduated,Natural	N.S.	1000	20000
YR1249	Pipette tip, 100 -1000 ul,Eppendorf with collar, Blue	N.S.	500	10000
YR1250	Pipette tip, 5 -200ul,MLA, Natural,1000/Zip bag,30000/Case,Non-sterile	N.S.	1000	20000
YR1251	Pipette tip, 100 -1000 ul,MLA,Natural,500/Zip bag,10000/Case,Non-sterile	N.S.	500	10000
YR1252	Pipette tip, 5 -200ul,Oxford, Natural,1000/Zip bag,30000/Case,Non-sterile	N.S.	1000	20000
YR1253	Pipette tip, 100 -1000 ul,Oxford,Blue,500/Zip bag,10000/Case,Non-sterile	N.S.	500	10000
YR1254	Pipette tip, 300ul,Extra Long 9cm,500/Zip bag,10000/Case,Non-sterile	N.S.	500	10000
YR1255	Pipette tip, 1000 ul,Extra Long 10.5cm,500/Zip bag,10000/Case,Non-sterile	N.S.	500	10000

Pipette Tips

Options available for Pipette Tips :

Variation

- (YR1251) Capacity: 100 -1000ul / Features: MLA, 500/Zip bag, 100 (-\$ 58.74),
- (YR1252) Capacity: 5 -200ul / Features: Oxford, 1000/Zip bag, (+\$ 38.08),
- (YR1241) Capacity: 5 -10ul / Features: Gilson crystal, graduated,
- (YR1253) Capacity: 100 -1000ul / Features: Oxford, 500/Zip bag, (+\$ 7.16),
- (YR1242) Capacity: 1 -200ul / Features: Gilson graduated, bevell,
- (YR1254) Capacity: 300ul / Features: Extra Long 9cm, 500/Zip bag (+\$ 583.92),
- (YR1243) Capacity: 1 -200ul / Features: Gilson graduated, bevell,
- (YR1255) Capacity: 1000ul / Features: Extra Long 10.5cm, 500/Zip (+\$ 735.84),
- (YR1244) Capacity: 5 -200ul / Features: Universal Gilson. Yellow,
- (YR1246) Capacity: 0.5-10ul / Features: Eppendorf crystal. Natural (+\$ 34.74),
- (YR1245) Capacity: 100 -1000ul / Features: Universal Gilson, gra (-\$ 33.81),
- (YR1247) Capacity: 5 -200ul / Features: Eppendorf with collar. Y,
- (YR1248) Capacity: 10 -300ul / Features: Universal Eppendorf gra (+\$ 8.10),
- (YR1249) Capacity: 100 -1000ul / Features: Eppendorf with collar (-\$ 121.45),
- (YR1250) Capacity: 5 -200ul / Features: MLA, 1000/Zip bag, 30000 (+\$ 8.10).