

Storage boxes for slide microscope



Model : YR1288

Storage boxes for slide microscope

Manufacturer : Kalstein

Price : \$ 441.69

Description :

Keeping everything organized; clean and secure is the key to run a reliable laboratory. One of the most used devices in that area is the microscope. This is a crucial factor to have reliable results when it comes to running tests. The slides you use to look at the samples have to be stored in a very safe place, so that they can stay away from humidity and more importantly, so that they do not break. Humidity can significantly affect the slides to look with the microscope. On the other hand, those slides need to be in a place where they do not receive any hits, and they need to be located where nothing break them.

These are high impact Acrylonitrile Butadiene Styrene (ABS) a common thermoplastic polymer, which has no true melting point, these means that it can be at an average temperature and it will not suffer any damage. First, this feature make this a very good material, because of its resistance, toughness, and specially heat resistance. Other advantages of this material is that it is also resistant to aqueous acids, alkalis, concentrated hydrochloric and phosphoric acids, alcohols and animal, vegetable and mineral oils. Besides the high impact ABS, this storage boxes have heavy walls, to keep humidity away and one of the most important features is that is unbreakable.

These boxes have a cork lining to protect them, which prevents any outer element to affect the inside, and comes in handy to prevent any damage to the stored slides. The measurements of each space for the slides are 25.4 X 76.2 mm (1" x 3").

Since the correct storage is one of the most important things in a laboratory, these boxes come with a numbered identification to correspond to the inventory sheet of the lab, which is located in the inside cover of the box for proper slide identification. Moreover, there is another feature, they come with five different coded colors, so that it is easier to identify among the other storage boxes. In addition, they are stackable, which makes them perfect for keeping them in shelves without losing the order.

We are offering you with five different capacity options for these storage boxes. The first one is a 12-places model with a box measurement of 80 mm x 54 mm x 33 mm (L x W x H). The second one is a 25-places model with a box measurement of 81 mm x 100 mm x 33 mm (L x W x H). The third model measures 195 mm x 95 mm x 33 mm (L x W x H) with a capacity for 50 places. The fourth one is a 100-places slip-on model with a box measurement of 190 mm x 170 mm x 33 mm (L x W x H). Finally, the last model measures 208 mm x 175 mm x 33 mm and it is a 100-place hinge pin-storage box.

Made from high impact ABS (with heavy walls), unbreakable and keep away

Storage boxes for slide microscope

from humidity.

The cork lining prevents damage to the stored slides (25.4 X76.2mm (1" x 3")).

Each slot is numbered to correspond to the slide inventory sheet on the inside cover of the box for proper slide identification.

Available in 5 colors coded boxes making the boxes easy to identify.

Stackable.

Parameters by Models

Storage boxes for slide microscope

Storage Boxes for Slide Microscope				
Model	Description	Sterility	Inner packing	Case Qty.
YR1288	Microscope slides box, 12place, 80 x 54 x 33mm (L x W x H)	N.S.	1/bag	240
YR1289	Microscope slides box, 25place, 81 x 100 x 33mm (L x W x H)	N.S.	1/bag	120
YR1290	Microscope slides box, 50place, 199 x 95 x 33mm (L x W x H)	N.S.	1/bag	120
YR1291	Microscope slides box, 100place, Slip-on, 199 x 170 x 33mm (L x W x H)	N.S.	1/bag	50
YR1292	Microscope slides box, 100place Hinge Pin, 208 x 175 x 33mm (L x W x H)	N.S.	1/bag	50

Storage boxes for slide microscope

Options available for Storage boxes for slide microscope :

Variation

- (YR1288) Dimensions: 80 x 54 x 33mm / Feature: 12 place,
- (YR1290) Dimensions: 199 x 95 x 33mm / Feature: 50 place (+\$ 63.10),
- (YR1291) Dimensions: 199 x 170 x 33mm / Features: 100 place. Sli (-\$ 96.98),
- (YR1292) Dimensions: 208 x 175 x 33mm / Features: 100 place. Hin (-\$ 91.14),
- (YR1289) Dimensions: 81 x 100 x 33mm / Feature: 25 place (-\$ 126.19).