

Portable Incubator



Model : YR0571

Portable Incubator

Manufacturer : Kalstein

Price : \$ 2,350.00

Description :

Having a portable incubator is one of the most comfortable advantages that a lab device can deliver. As you may know, incubators are an essential device in a laboratory when it comes to cellular biology experiments, microbiology, and cellular cultures. These are very important devices because you need them to maintain the right temperature, humidity and other factors involved in experiments that imitate the living process. Kalstein products provides you the ideal constant temperature and humidity incubator so that you can rely on the storage of any kind of culture you want to preserve.

Portable incubator

A small and lightweight unit is the ideal device you need to have when it comes to personal uses. For example, if you need to take lab-like conditions to the field on an onsite testing of bacterial contamination, you must have a portable incubator. They are a necessary item to create reliable and repeatable conditions in the field. These models are trust-worthy devices with stable temperatures that guarantee a stable environment for the incubation of bacteriological samples.

These units are ideal for accurate incubation of total coliform, fecal coliform and E. coli testing. In addition, you can use them to incubate samples at temperatures ranging from 5°C to 60 °C, when it comes to performing experiments at different setting temperatures.

Constant temperature incubators

These models have terrific features that guarantees the safety and quality of the cultures you want to maintain. These models have a stainless steel inner chamber, making them high quality devices for your laboratory, being this one of the most desired features in this kind of devices. Another great feature that assures the complete sealing of the incubators' door is the silicon magnetic sealing, providing reliable results for your laboratory. In addition, you should also know that these models have a PID control, which means that this controller automatically applies accurate and responsive correction to a control function. Moreover, another important safety feature is the over temperature alarm.

Since you are purchasing high quality products for your laboratory, it is important that you know the models from the series. There is the YR0571 and the YR0572 model; their measurements are 200 mm x 250 mm x 200 cm for the internal size and 340 mm x 260 mm x 360 mm for the external

Portable Incubator

size. Moreover, their power consumption is 100W and their weight is 9.7 Kg.

Features of the Portable Incubator

- Stainless steel inner chamber
- Silicon magnetic sealing
- PID control
- Over temperature alarm

Parameters by Model

Portable Incubator

Floor-standing High-Speed Refrigerated Centrifuge

Model:	YR0571	YR0572
Power supply:	AC220V 50Hz	DC12V/AC220V
Temperature range (°C):	RT +5~45°C	
Consumption:	100W	
Internal size (mm):	200*250*200	
External size (mm):	340*260*360	
Packing size (mm):	400*330*425	
Gross weight(kg):	9.7	

Portable Incubator

Options available for Portable Incubator :

Variation

(YR0571) 220VAC,

(YR0572) 220VAC/12VDC (+\$ 229.61).