

## Transfer Pipette 4ml Capacity - Short Stem



**Model : YR1143**

**Transfer Pipette 4ml Capacity - Short Stem**

**Manufacturer : Kalstein**

**Price : \$ 336.53**

**Description :**

Many items are an important part of an inventory in a lab. Without them, testing samples and working in any medical environment would not be impossible. Kalstein products is offering you the possibility of having the high-end tools and products for your laboratory, so that you can rely and trust that the work you are doing, is the best.

One of the most common tools and devices inside any lab is the transfer pipette. This is a lab tool, commonly used in chemistry, biology and medicine to transport a measured volume of liquid, often as a media dispenser. There are several designs and purposes with different levels of accuracy and precision. In this case, we are presenting you with the 4 ml capacity transfer pipette with short stem.

### Transfer pipette 4 ml

This kind is the most basic pipettes in the area. You can use it for rough measurements only; however, it is still very important that you follow a standard pipetting technique. This kind of pipette is a single piece of plastic pipette with a short stem, which means that you will have to be more careful when you handle any kind of sample.

### Transfer pipette function

About the functions, it is recommended that the user aspirates the liquid at a 90 degree angle, dispense at a 45 degree angle and touch off to make sure all the liquid is dispensed. These recommendations are a very important advice to follow to ensure the precision in your sample testing; this way guarantees less mistakes in the results, making them more reliable. The uses ranges from general purpose to sensitive lab research.

The particular pipette we are presenting you has a 4ml capacity. The total length of the transfer pipette is 83 mm. The stem diameter is 2.5 mm. In addition, it has a standard bulb, which means that the draw is 3.5 ml. Finally, the dropper/ml is 22.

Finally, there is an important thing to know about the packaging, the transfer pipettes can come as 500 unit per box and 3000 per case, this category is

## Transfer Pipette 4ml Capacity - Short Stem

the non-sterile type. There is also the option of packaging 1 per peel pack, 500 per case, and 5 per peel pack and 500 per case. These last two options have gone through the Ethylene Oxide Sterilization process.

- Capacity:4ml
- Total Length: 83mm
- Dropper/ml: 22
- Stem Diameter: 2.5mm
- Bulb Draw:3.5ml

Parameters by Model

## Transfer Pipette 4ml Capacity - Short Stem

**Transfer pipette, 4ml Capacity - Short Stem**

Model	Description	Sterility	Case Qty.
YR1143	500/Box,3000/Case,Nonsterile	Non-	3000
YR1144	1/Peel Pack,500/Case,Sterile	EO	500
YR1145	5/Peel Pack,500/Case,Sterile	EO	500

## Transfer Pipette 4ml Capacity - Short Stem

### Options available for Transfer Pipette 4ml Capacity - Short Stem :

#### Variation

- (YR1143) 500/Box 3000/Case Nonsterile,
- (YR1144) 1/Peel Pack 500/Case Sterile (-\$ 141.39),
- (YR1145) 5/Peel Pack 500/Case Sterile (-\$ 149.57).