

Petri Dish



Model : YR1214

Petri Dish

Manufacturer : Kalstein

Price : \$ 163.59

Description :

When it comes to the study of microorganisms, it is very important that the object of study is isolated from other species. That is why having the right laboratory lab ware is a very important aspect in order to achieve reliable results. A petri dish or cell-culture dish is a shallow cylindrical round lidded dish used to culture cells. Those can be bacteria or any kind of viruses.

Petri dish

This petri dish model is the ideal tool you must have in your laboratory in order to culture microorganisms on media. It is a virgin polystyrene petri dish, which means that they are manufactured with a material that is suitable for all kinds of thermoformed packaging, providing the necessary safety for any containers. This is a very important feature because it assures confidence to every dish you will use in your daily workflow.

We know that quality is a very important aspect to you. We also care about that, and that is why this petri dish model has 100,000-grade clean room guarantee quality with the Management System EN ISO13485. In addition, these petri dishes are CE marked in accordance with the European Directive 98/79/EC.

It is important for you to know that they are for in vitro use only, which means that they are strictly to perform tests with microorganisms, cells, or biological molecules outside their normal biological context. There is another physical aspect you should know about this petri dish, and that is that it has mirror-finishing that ensures high clarity, being this a very important feature when it comes to working on a daily basis in a microbiology lab. Moreover, you should know that the 90 mm dish model is designed for an automatic plate pourer.

There are several versions or models of this petri dish. The first one measures 150 mm x 15 mm and it has triple vents. The second model measures 130 mm x 130 mm and it is square-shaped with triple vents. In addition, there is three petri dish models that measure 90 mm x 15 mm, one of them only has triple vents, the other one has 2 compartments and the third one has 3 compartments. The last two also have triple vent. Moreover, one of the versions measure 65 mm x 15 mm and it has counting areas with triple vents. In addition, there is another version that measures 55 mm x 15 mm with triple vents and finally, the last one measures 35x15mm and it also has triple vents.

- Ideal for culture micro-organisms on media.
- Manufactured from virgin Polystyrene assuring confidence to every dish.
- Management System EN ISO13485 with SOPs and 100,000grade

Petri Dish

clean-room guarantee quality.

- CE marked in accordance with European Directive 98/79/EC
- Invitro use only.
- Mirror-finishing ensures high clarity.
- 90mm Dish designed for automatic plate pourers.

Parameters by Model

Petri Dish

Petri Dish

Model	Description	Sterility	Inner packing	Case Qty.
YR1214	Petri dish, 150x15mm, triple vents.	EO	1/bag	100
YR1215	Petri dish, 130x130mm, square, triple vents	EO	1/bag	150
YR1217	Petri dish, 90x15mm, 2 compartments, triple vents	EO	10/bag	500
YR1218	Petri dish, 90x15mm, 3 compartments, triple vents	EO	10/bag	500
YR1219	Petri dish, 65x15mm, counting areas, triple vents	EO	10/bag	1000
YR1220	Petri dish, 55x15mm, triple vents.	EO	10/bag	1000
YR1222	Petri dish, 35x15mm, triple vents.	EO	10/bag	3000

Petri Dish

Options available for Petri Dish :

Variation

- (YR1219) Petri dish 65x15mm Counting areas | Triple vents (+\$ 116.00),
- (YR1220) Petri dish 55x15mm Triple vents (+\$ 23.37),
- (YR1214) Petri dish 150x15mm Triple vents,
- (YR1215) Petri dish 130x130mm Square | Triple vents (+\$ 81.80),
- (YR1217) Petri dish 90x15mm 2 compartments | Triple vents (-\$ 2.92),
- (YR1218) Petri dish 90x15mm 3 compartments | Triple vents (+\$ 93.00).