

Refractometer



Model : YR01768

Refractometer

Manufacturer : Kalstein

Price : \$ 92.87

Description :

Refractometer is precision optical instrument that operates on the principle- as the concentration or density of a solution increases, its refractive index changes proportionately. The refractive angle measured by refractometer registers on the scale. The larger the concentration of solute in the solution, the higher the reading on the scale.

Parameters by Models

Refractometer

Refractometer

Model	Descripción
YR01768	With ATC for low concentration 0-10% Brix div 0,1%
YR01769	With ATC for low concentration 0-20% Brix div. 0,1%

Refractometer

Options available for Refractometer :

Variation

(YR01768) With ATC for low concentration 0-10% Brix div 0,1%,
(YR01769) With ATC for low concentration 0-20% Brix div. 0,1%.