

Refractometer



Model : YR01770
Refractometer
Manufacturer : Kalstein

Price : \$ 92.87

Description :

Refractometer is precision optical instrument that operates on the principle- as the concentration or density of a solution increases, its refractive index changes proportionately. The refractive angle measured by refractometer registers on the scale. The larger the concentration of solute in the solution, the higher the reading on the scale

Parameters by Models

Refractometer

Réfractomètre

Model	Description
YR01770	Avec ATC 0-32% Brix div. 0,2%
YR01771	avec ATC 28-62% Brix div. 0,2%
YR01772	avec ATC 40-82% Brix div. 0,5%

Refractometer

Options available for Refractometer :

Variation

(YR01770) With ATC 0-32% Brix div. 0,2%,
(YR01771) with ATC 28-62% Brix div. 0,2%,
(YR01772) with ATC 40-82% Brix div. 0,5%.