

## Refractometer



**Model : YR01773**  
**Refractometer**  
**Manufacturer : Kalstein**

**Price : \$ 92.87**

**Description :**

Refractometer is precision optical instrument that operates on the principle- as the concentration or density of a solution increases, its refractive index changes proportionately. The refractive angle measured by refractometer registers on the scale. The larger the concentration of solute in the solution, the higher the reading on the scale

Parameters by Models

## Refractometer

### Refractometer

Model	Descripción
YR01773	with ATC for honey 58-92% Brix div.0.5%, 12-27% water div.0,5%
YR01779	with ATC for sodium chloride 0-28% Brix div. 0,2%

## Refractometer

### Options available for Refractometer :

#### Variation

(YR01773) with ATC for honey 58-92% Brix div.0.5%, 12-27% water div.0,5%,  
(YR01779) with ATC for sodium chloride 0-28% Brix div. 0,2%.