

Refractometer



Model : YR01774
Refractometer
Manufacturer : Kalstein

Price : \$ 92.87

Description :

Refractometer is precision optical instrument that operates on the principle- as the concentration or density of a solution increases, its refractive index changes proportionately. The refractive angle measured by refractometer registers on the scale. The larger the concentration of solute in the solution, the higher the reading on the scale.

Parameters by Models

Refractometer

Refractometer

Model	Descripción
YR01774	with ATC 0-50% Brix div. 0,5%
YR01775	with ATC 0-80% Brix div. 0,5%
YR01776	with ATC 0-90% Brix div. 0,5%

Refractometer

Options available for Refractometer :

Variation

(YR01774) with ATC 0-50% Brix div. 0,5%,
(YR01775) with ATC 0-80% Brix div. 0,5%,
(YR01776) with ATC 0-90% Brix div. 0,5%.