

## Refractometer



**Model : YR01780**

**Refractometer**

**Manufacturer : Kalstein**

**Price : \$ 92.87**

**Description :**

Refractometer is precision optical instrument that operates on the principle- as the concentration or density of a solution increases, its refractive index changes proportionately. The refractive angle measured by refractometer registers on the scale. The larger the concentration of solute in the solution, the higher the reading on the scale

Parameters by Models

## Refractometer

### Refractometer

Model	Descripción
YR01780	with ATC 0-10% Brix div.0.1%, 0-100‰ salinity div.1%
YR01781	with ATC 0-32% Brix div.0.2%, 0-28‰ salinity div.0,2%

## Refractometer

### Options available for Refractometer :

#### Variation

(YR01780) with ATC 0-10% Brix div.0.1%, 0-100% sanility div.1%,  
(YR01781) with ATC 0-32% Brix div.0.2%, 0-28% sanility div.0,2%.