

## Refractometer



**Model : YR01782**  
**Refractometer**  
**Manufacturer : Kalstein**

**Price : \$ 92.87**

**Description :**

Refractometer is precision optical instrument that operates on the principle- as the concentration or density of a solution increases, its refractive index changes proportionately. The refractive angle measured by refractometer registers on the scale. The larger the concentration of solute in the solution, the higher the reading on the scale

Parameters by Models

## Refractometer

### Refractometer

Model	Descripción
YR01782	with ATC for protein 0-12g/dl div.0-0.2%/dl, 1.000-1.040sg div.0.002sg
YR01783	with ATC for protein 0-12g/dl div.0-0.2%/dl, 1.000-1.050sg div.0.002sg,1.3330-1.3600RI div.0.0005RI

## Refractometer

### Options available for Refractometer :

#### Variation

(YR01782) with ATC for protein 0-12g/dl div.0-0.2%g/dl, 1.000-1.040sg div.0.002sg,  
(YR01783) with ATC for protein 0-12g/dl div.0-0.2%g/dl, 1.000-1.050sg div.0.002sg, 1.3330-1.3600RI  
div.0.0003RI.