

SafeBox mini range



Model : 00250-01-08-00

SafeBox mini range

Manufacturer : SafeBox

Price : \$ 0.00

Description :

SafeBox mini range

The Mini range is a set of tamper-proof pocket collectors used for transporting and disposing of needles, syringes and other sharp objects. The ergonomic design gives you an easy and safe work practice. They are known for their flexibility and are in sizes of 0.25, 0.5 and 1 liter. They are small enough to be carried with you on your daily tour. The Mini range is designed with extra wall thickness to prevent punctures. All sizes are stackable both empty and full. It is perfect for diabetics because it is equipped with an insulin pen selector.

In the lid there is also a QuickFit alternative and a universal disconnect switch. The container can be locked temporarily or permanently. The Mini range has a visible fill line that indicates when the container is full. All sizes are fully autoclavable.

Advantages

- UN and ADR standards are validated
- stackable
- Marking for maximum fill level
- autoclavable
- Temporary closure
- Made of ecological polypropylene
- Hypodermic needle puller in the lid
- Designed with different types of support: wall, flat surface and trolley
- Extra wall thickness to prevent punctures

*PRESENTATION 1000 PIECES

SafeBox mini range

REFERENCES, CAPACITIES, DIMENSIONS AND WEIGHTS

Model:	REF: 00250-01-08-00 Mini de 0,25 liter	REF: 00500-01-08-00 Mini de 0,5 liter	REF: 01000-01-08-00 Mini de 1 liter
Weight (g):	34,5	59	94,5
Nominal Capacity (liter):	0,24	0,5	1
Maximum capacity 95% of overflow (liter):	0,22	0,46	0,95
Useful capacity (liter):	0,2	0,4	0,8
Dimension H x D (NET mm) :	124x60x34	144,6x81,5	256x81,5
Minimum thickness of walls (mm):	1,47	1,81	1,55
Total height with cover (mm):	1,3	166	276

SafeBox mini range

Options available for SafeBox mini range :

Capacity (Liters)

Mini 0,25 liter,
Mini 0,25 liter,
Mini 0,5 liter (-\$ 30.80),
Mini 1 liter (+\$ 33.60).

Labeled

Mini 1 L Labeled (+\$ 77.00),
Mini 0.25 L Labeled (+\$ 44.80),
Mini 0.5 L Labeled (+\$ 14.00).