

EDTA K2 Tube



Model : VG30031S

EDTA K2 Tube

Manufacturer : -

Price : \$ 0.00

Description :

It is coated with EDTAK2. It is widely used in clinical haematology and various kinds of blood cell test instrument. Used for direct sampling and analyzers without actually being opened.

EDTA K2 Tube

Capacities and dimensions

Product Name	EDTAK2 Tube	
Material	Pet, Glass	
Draw Volume	Pet: 13x75MM:1ML;2ML;3ML;4ML 13x100MM: 4ML;5ML;6ML 16X100MM:7ML;8ML;9ML	Glass: 13x75MM:1ML;2ML;3ML;4ML;5ML 13x100MM: 4ML;5ML;6ML;7ML 16X100MM:8ML;9ML;9.5ML;10ML 16X100MM: 6ML;7ML ;8ML;8.5ML
Cap Color	LAVENDER	
Closure Type	SAFETY	
Package	100 PCS/Tray; 1200 PCS/Carton	
Measurement	Carton Measurement of 13x75MM: 515x345x205MM/1200PCS Carton Measurement of 13x100MM: 515x345x255MM/1200PCS Carton Measurement of 16x100MM: 590x392x262MM/1200PCS	
Shelf Life	1 Year/PET, 2 Years/Glass	
Others	Storage in room temperature. Single use only! Dispose it once used!	

EDTA K2 Tube

Options available for EDTA K2 Tube :

Capacity (mm) Pet

13x100mm,
13x75mm,
16x100mm.

Capacity (mm) Glass

13x100mm,
13x75mm,
16x100mm.