

SafeBox Prime Range



Model : 04000-01-08-00

SafeBox Prime Range

Manufacturer : -

Price : \$ 0.00

Description :

The Prime range is a very flexible and innovative set of high quality square containers. They are specifically designed for a secure work environment in a hospital or healthcare facility. The Prime range is designed with extra wall thickness to prevent punctures. All sizes are stackable both empty and full. The manifolds are all equipped with a universal needle switch, universal blade disconnect, QuickFit alternative, Luer lock switch and slide switch. They can be locked temporarily or permanently. The permanent closure is indicated by a latch drawn on the flap. The Prime line has a visible fill line that indicates when the container is full. All sizes are fully autoclavable and equipped with handle for easy transport. Available sizes are 4 liters, 6 liters and 11 liters.

ADVANTAGES

- UN and ADR standards are validated
- stackable
- Marking for maximum fill level
- autoclavable
- Temporary closure
- Made of ecological polypropylene
- Hypodermic needle puller in the lid
- Designed with different types of support: wall, flat surface and trolley
- Extra wall thickness to prevent punctures

MAIN FEATURES

-
- Available sizes: 4 -6 and 11 liters
- Standard color: Yellow
- Material: PP copolymer

***PRESENTATION 1000 PIECES**

SafeBox Prime Range

Capacities, dimensions and weights

Model:	04000-01-08-00	06000-01-08-00	11000-01-08-00
	Prime 4 liters	Prime 6 liters	Prime 11 liters
Weight (g):	232,5 g	324 g	531,5 g
Nominal capacity (liter):	4 litres	6 litres	11 litres
Maximum capacity at 95% filling:	3,4 litres	5,1 litres	9,25 litres
Minimum thickness of the walls:	3,2 mm	4,8 mm	8,8 mm
Dimension H x D (NET mm):	138x280x190 mm	202x280x190 mm	323x280x190 mm
Weight of the lid:	127	127	127
Cover size H x D:	33x258x305	33x258x305	33x258x305
Total height with lid:	145 mm	210 mm	330 mm

SafeBox Prime Range

Options available for SafeBox Prime Range :

Capacity (Liters)

Guardian 4 liters,
Guardian 6 liters (+\$ 462.00),
Guardian 11 liters (+\$ 924.00).

Labeled

Guardian 4L Labeled (+\$ 75.60),
Guardian 6 L Labeled (+\$ 536.20),
Guardian 11 L Labeled (+\$ 999.60).