

Benchtop DO/BOD/OUR/SOUR Meter

Benchtop DO/BOD/OUR/SOUR Meter



Model : YR01831

Benchtop DO/BOD/OUR/SOUR Meter

Manufacturer : Kalstein

Price : \$ 2,075.00

Description :

YR01831 Benchtop DO/BOD/OUR/SOUR Meter

Professional benchtop dissolved oxygen meter, selectable 3 measurement modes. The meter is suitable for measuring the dissolved oxygen, BOD, Oxygen Uptake Rate and Specific Oxygen Uptake Rate of the samples.

Features:

- Professional benchtop dissolved oxygen meter is equipped with a 7 inch TFT display.
- 1 or 2 points push-button calibration using the air-saturated water or zero oxygen solution.
- Salinity and barometric pressure compensations ensure the accurate DO measurements.
- Selectable testing time, beginning/ending DO are used for the OUR and SOUR calculations.
- Limit alarm alerts when reading exceeds range.
- Calibration Due Alarm prompts user to calibrate the meter regularly.
- Stability indicator shows when a measurement is recognized as stable.
- Auto-Read feature senses and locks the measurement endpoint.
- Interval Readings automatically send the measured data to the computer or printer.
- Password protection prevents the unauthorized calibration and settings.
- Expanded memory stores and recalls up to 1000 data sets.
- Stored readings can be transferred into the computer by USB communication interface.
- Reset feature automatically resumes all settings back to factory default options.

[Specifications](#)[Order Information](#)[Optional Accessories](#)[Packaging](#)