

## Mini Vortex Mixer



**Model : YR01887**

**Mini Vortex Mixer**

**Manufacturer : Kalstein**

**Price : \$ 200.00**

**Description :**

Simple device used commonly in laboratories to mix small vials of liquid.

They may be used to mix the reagents of an assay or to mix an experimental sample and a diluent. perfectly smooth protein chemical or sensitive culture applications;

In cell culture and microbiology laboratories they may be used to suspend cells. In a biochemical or analytical laboratory they may be used to mix reagents of an assay or to mix an experimental sample and dilutant

Product Name: Mini Vortex Mixer.

Rated frequency [Hz]: 50-60

Fixed speed: 3000rpm

Application: It is the necessary equipment when the laboratory staff mixing for a variety of liquid. This is an ideal fixture for small sample test, convenient to use in many environment (indoor/outdoor). In cellculture and microbiology laboratories they may be used to suspend cells. In a biochemical or analytical laboratory.

Rated voltage [ V AC]: 100-240

Other Configuration:

- Orbital diameter 5.5 [mm]
- tubes size diameter 30mm can be placed.
- Touch function.
- Open Lid Automatically.