

---

*Laboratory Lines > Dry Bath Incubators > Mini Dry Bath Incubator With heating lid without block included*

## **Mini Dry Bath Incubator With heating lid without block included**

## Mini Dry Bath Incubator With heating lid without block included



**Model : YR01882**

**Mini Dry Bath Incubator With heating lid without block included**

**Manufacturer : Kalstein**

**Price : \$ 565.00**

**Description :**

Mini Dry Bath Incubator With heating lid without block included.

YR01882 Mini dry bath incubator oem products dental apparatus

Mini Dry Bath Incubator fits conveniently on every bench with its very small footprint. The digital microprocessor ensures temperature accuracy of the metal block and quick heating saves your valuable time. This small but powerful unit is suitable for a range of applications: Restriction digests, denaturing DNA, NUB, melting agar, coagulation studies, in-situ hybridisation, and Hot Start thermocycler reactions. The unit is supplied with a transparent lid for temperature uniformity across all tubes.

**Feature:**

- LED display for easy reading
- Heating lid
- Temperature is up to 100°C
- Simultaneous display of set/actual time and temperature
- 3D digital wheel adjustment of speed and temperature
- High temperature precision is up to  $\pm 0.5^{\circ}\text{C}$
- Programmable operation can store up to 5 programs
- Low operation voltage (12v) can be used directly with car charger
- Various applications with different blocks
- Block for tubes: 0.2ml, 0.5ml, 2ml, 15ml, 50ml tube

**Parameters by Model**

## Mini Dry Bath Incubator With heating lid without block included

<b>Product Name</b>	<b>YR01882</b>
<b>Temperature Control Range</b>	RT + 5°C~100°C
<b>Timing Range</b>	0~99h59min
<b>Temp. Control Accuracy</b>	±0.5°C
<b>Display Accuracy</b>	±0.1°C
<b>Heating Time (25°C to Max. Temp.)</b>	≤15min
<b>Dimension(L x W x H) mm</b>	164 x 116 x 98
<b>Net Weight</b>	0.90kg
<b>Rated voltage [ V AC]</b>	12V