

Binocular Microscope



Model : YR-0237

Binocular Microscope

Manufacturer : Kalstein

Price : \$ 1,975.00

Description :

Biological Microscope

Description

Excellent image quality with infinitive optical system.

Comfortable operating with ergonomic structure.

Better illumination with sliding-in centerable condenser, Convenient replacement with phase contrast and dark field condenser.

Features?

- Extra Wide Field Eyepiece EW10X/20 with Diopter Adjustment as Optional.
- Comfortable Observation with Eyepoint and interpupillary adjustable.
- Sliding-in Centerable Condense Adjustable.
- Handle for easy carrying.
- Lamp Replacement by Direct Removing Housing.
-

Binocular Microscope

Parameters by Models

Binocular Microscope

Model	YR0237		
Viewing Head	Compensation Free binocular head, 30° inclined, Interpupillary 48-75mm	●	
	Compensation Free trinocular head, 30° inclined, Interpupillary 48-75mm	○	
Eyepiece	Wide Field Eyepiece WF 10×/18	○	
	Extra Wide Field Eyepiece EW10×/20 with Diopter Adjustment	●	
Objective	Infinitive Semi-plan Achromatic Objectives	4x	●
		10x	○
		40x	●
		100x	●
Nosepiece	Quadruple Nosepiece	●	
Stage	Double Layers Mechanical Stage 140mm×140mm/ 75mm×50mm	●	
Condenser	Sliding-in centerable condenser NA12.5	○	
Focusing	Coaxial Coarse & Fine Adjustment, Moving Range 20mm	○	
Illumination	6V/ 20W Halogen Lamp, Brightness Adjustable	○	
Phase Contrast Kit		○	
Dark Field Attachment		○	

Note: ● Standard Outfit, ○ Optional

Binocular Microscope

Options available for Binocular Microscope :

Variation

(YR0237)BINOCULAR,
(YR0237-1)TRINOCULAR (+\$ 480.00).