

## Microplate Reader and Washer STD Config



**Model : YR02208**

**Microplate Reader and Washer STD Config**

**Manufacturer : Kalstein**

**Price : \$ 6,000.00**

**Description :**

**MICROPLATE READER & WASHER STD CONFIG**

**Features:**

- Washes flat, U and V-bottom strips and plates.
- Small residual volume, multiple aspiration ways.
- Multi-channel, incubator optional.
- Liquid level detection and alert function.

**Parameters by Model**

## Microplate Reader and Washer STD Config

<b>Model</b>	<b>YR02208</b>
<b>Plate Mode</b>	96- or 48- well plate or strip
<b>Residual volume</b>	<3 ul
<b>Washing volume</b>	10-3000 ul/well, 1 ul by step
<b>Washing cycle</b>	Maximum 99 cycles
<b>Strip washing</b>	1 to 12 strips programmable
<b>Soaking time</b>	0-24 hrs
<b>Shaking time</b>	0-24 hrs
<b>Display</b>	Large LCD
<b>Wash Channels</b>	5 channels, 1 or up to 3 for wash
	1 forrinse and 1 for waste
<b>Memory</b>	99 washing protocols
<b>Interface</b>	RS-232
<b>Net Weight</b>	AC 110V/220V, 50Hz
<b>Dimensions LxWxH (mm)</b>	450*390*190
<b>Incubator temperature (optional)</b>	25 ,30 ,37

## Microplate Reader and Washer STD Config

### Options available for Microplate Reader and Washer STD Config :

#### Variation

(YR02086) CONFIG + INCUBATOR (+\$ 750.00).

#### Variation

(YR02208) MICROPLATE READER & WASHER STD CONFIG,

(YR02209) MICROPLATE READER & WASHER STD CONFIG + INCUBATOR (+\$ 750.00).