

## Binocular Microscope



**Model : YR0246**  
**Binocular Microscope**  
**Manufacturer : Kalstein**  
**Price : \$ 4,725.00**  
**Description :**

Series Stereomicroscope

Description:

Excellent Image Quality & Optical Performance with Infinity Parallel Optical System.

Comfortable Operation with Ergonomic Design Principle.

With LED light for both incident and transmitted illumination, Providing even illumination and life expectancy can reach 6000 hours.

A Wide Range of Accessories to be Equipped with the Microscope for Multi-purposelt is an Outstanding Assistant in Industrial, Medical and Scientific Research Areas.

Parameters by Models

## Binocular Microscope

Model		YR0246	YR0247	YR0248
Drawing Tube	Binocular Head, 20°Inclination	•	•	•
	Tilting Binocular Head, 0°-35° Inclination			
	Photo Attachment (one photo tube)			
	Photo Attachment(two photo tube)	○	○	○
Eyepiece	EW10×/22	•	•	•
	WF15×/16	○	○	○
	WF20×/12	○	○	○
Zoom Objective	0.8×-5×	•		
	0.8×-6.4×		•	
	0.8×-8×			•
Objective	Plan Achromatic Objective 1×	•	•	•
	Plan Apochromatic Objective 0.5×	○	○	○
	Plan Achromatic Objective 2×	○	○	○
	Achromatic Objective 0.3×			
Zoom Ratio	1: 6	•		
	1: 8		•	
	01:10			•
Working Distance	78mm	•	•	•
Focusing	Focusing Range 105mm	•	•	•
Illumination	Transmission or Reflection LED Illumination, Brightness Adjustable	•	•	•
	Coaxial Illumination	○	○	○

Note: • Standard Outfit, ○ Optional

## Binocular Microscope

### Options available for Binocular Microscope :

#### Variation

(YR0246),  
(YR0247) (+\$ 475.00),  
(YR0248) (+\$ 2,835.00).