

Coagulation Analyzer with electrical pipette



Model : YR02218

Coagulation Analyzer with electrical pipette

Manufacturer : Kalstein

Price : \$ 6,000.00

Description :

COAGULATION ANALYZER

With electrical pipette

Features:

- 4 channels (4 durable LED detectors for 4 different parameters analysis).
- QC and calibration program incorporated.
- Large memory, up to 10,000 results can be stored.
- Scattered light principle and percentage analysis.
- Easy Windows operation system with touch screen or mouse, large LCD display.
- Advanced scattered light principle and percentage analysis ensure accurate results.
- Low reagent consumption.
- Multiform result output including patient comprehensive report
- Auto power-on self-check.

Parameters by Model

Coagulation Analyzer with electrical pipette

Model	YR02218
Principle	Scattered light with percentage analysis
Test Channel	4 Channels
Light Source	Durable LED detector
Test Time	Normally 20-60 secs, Maximum time up to 600 secs
Sample Position	24 Positions
Reagent Position	6 Positions (1 magnetic mixer included)
Minimum Reagent Consumption	25 ul for TT, FIB, PT, APTT
Memory	10000 Test results
Interface	RS-232, PSTN (Optional)
Display	6"LCD
Input	Touch panel or mouse
Output	External printer (optional)
Net Weight	8Kg
Dimensions LxWxH (mm)	410*310*160
Power Supply	AC 110V-220V ±10%, 50/60Hz