

Patient Monitor New Version



Model : YR02178

Patient Monitor New Version

Manufacturer : Kalstein

Price : \$ 3,200.00

Description :

PATIENT MONITOR

10.4" screen, ECG, SpO2, RESP, NIBP, TEMP & PR AND 12" screen, ECG, SpO2, RESP, NIBP, 2*TEMP & PR

Multi-parameter Patient Monitor

Good value for money; 10.4 inch screen for YR02179, 12.1 inch screen for YR02178; multiple parameters available in a compact body; no-fan design.

Features

- Adult, kid and newborn monitoring are supported;
- Slim body with a hidden handle, easy to move;
- Visual and audio alarm with triple colours alarm light;
- Lithium battery inside;
- No-fan design, match the requirements of multi-patient ward;
- I Able to save data when power supply fails.

Parameters by Model

Patient Monitor New Version

Model	YR02178
NIBP	
Method	Oscillometric
Measurement modes	Manual/Automatic/Continuous
Displayed Parameters	Systolic, Diastolic, and MAP (mean arterial pressure)
Range of Systolic pressure	Adult: 40-270mmHg
	Pediatric: 40-200mmHg
	Neonate: 40-135mmHg
Range of Diastolic pressure	Adult: 10-210mmHg
	Pediatric: 10-150mmHg
	Neonate: 10-100mmHg
Range of Mean pressure	Adult: 20-230mmHg
	Pediatric: 20-165mmHg
	Neonate: 20-110mmHg
Accuracy	±6mmHg
Resolution	1mmHg
Over pressure protection	Double safety protection
SPO2	
Measurement range	0%-100%
Resolution	1bpm
Accuracy	±2%(70%-100%);unspecified(0-69%)
Alarm range	0%-100%
Pulse rate range	20 - 254 bpm
Resolution	1 bpm
Accuracy	±3 bpm
Alarm range	0 - 254 bpm
Temperature	
Measurement range	0-50°C
Resolution	0.1°C
Accuracy	±0.1°C
Channel	T1, T2, TD
Respiration	
Method	Thoracic impedance
Measurement range	Neonate/Pediatric: 0-150RPM
	Adult: 0-120RPM
Apnea alarm	No alarm,10s,15s,20s,25s,30,35s,40s, Selectable
Alarm	Audible and visual alarm, alarm events recallable
General Specifications	
Display	12.1"Color TFT LCD with LED background light, 800x600 resolution OxyCRG dynamic view display
Trace	7 traces

Patient Monitor New Version

Options available for Patient Monitor New Version :

Variation

(YR02178) 10.4" screen, ECG, SpO2, RESP, NIBP, TEMP & PR,

(YR02179) 12" screen, ECG, SpO2, RESP, NIBP, 2*TEMP & PR (-\$ 100.00).