

Pulse Oximeter



Model : YR02254
Pulse Oximeter
Manufacturer : Kalstein
Price : \$ 0.00

Description :

Introduction

The YR02254 Pulse Oximeter adopts Photoelectric Oxyhemoglobin Inspection Technology in accordance with Capacity Pulse Scanning & Recording Technology, which can be used to measure human oxygen saturation and pulse rate through finger. The device is suitable for being used in family, hospital, oxygen bar, community healthcare and physical care in sports, etc. (It can be used before or after exercise, but not recommended to use during exercising).

Features

- 1.3"Color OLED Display.
- Menu Operation.
- Battery Status Indication.
- Rechargeable Lithium Battery.
- Audible and Visible Alarm For Over-limit and Low Battery.
- Continuous 24h Storage for SPO2 and PR Value, Data Stored and Real-time data can be uploaded to PC