

Visual Electronic Stethoscope



Model : YR02266

Visual Electronic Stethoscope

Manufacturer : Kalstein

Price : \$ 300.00

Description :

Introduction

Visual Electronic Stethoscope is a kind of multi-functional electronic stethoscope. It takes the latest sensor technology to monitor human's Heart Sound and Lung Sound. It can detect real-time ECG waveform and calculate the heart rate, the pulse oxygen Saturation and pulse rate can also be monitored through SpO2 probe. The product is applicable for family, hospital, private Clinique, community medical treatment.etc.

Features

- 1.8"LCD.
- Display of ECG Waveform, HR, SPO2, PR, Pulse Waveform.
- Earphone Output and Adjustable 16-step Volume.
- Alarm for SPO2 and PR Exceeding Limits.
- Low Power Consumption.