

Ultrasound Scanner



Model : YR02267
Ultrasound Scanner
Manufacturer : Kalstein

Price : \$ 7,500.00

Description :

Introduction

The notebook type YR02267 is a full digital B-Ultrasound diagnostic system. It adopts embedded operating system, which greatly optimizes the product performance. The system is taken more conveniently for its high effective data processing ability, pop-up menu and keyboard design. It includes rich measuring software packages that satisfy the clinic diagnostic need fully.

Features

- 10"LCD
- Laptop Device with Embedded Operating System.
- Can be applicable for Examination of Abdomen, Gynecology and Obstetrics, Cardiology, Small Organs and Urology.
- Adopt THI Technology, reports can be printed and exported.
- Images can be stored in the device or U Disk.
- Be Compatible with Ink Jet Printer, Laser Printer and Video Printer.