

## Urine Analyzer



**Model : YR02272**  
**Urine Analyzer**  
**Manufacturer : Kalstein**  
**Price : \$ 0.00**  
**Description :**

### Introduction

YR02272 Urine Analyzer is a high-precision, intellectual instrument which is researched and developed basing on modern optics, electronics, computer science and other advanced technologies for clinical inspection of urine. GLU, BIL, SG, KET, BLD, PRO, URO, NIT, LEU, VC and PH in urine can be tested by using it with special test strips. It can be widely used in various medical and health departments as one of the main clinical laboratory instruments.

### Features

2.8"Color LCD.

Built-in Thermal Printer.

Standard RS232 and USB Interface.

Power Supply: AC.

Memory 1000 Sample Data.

Optional Units: International Unit, Conventional Unit and Symbol System.

Test Items: GLU, BIL, SG, KET, BLD, PRO, URO, NIT, LEU, VC, PH.