

SpO2 Simulator



Model : YR02281

SpO2 Simulator

Manufacturer : Kalstein

Price : \$ 0.00

Description :

Introduction

YR02281 SpO2 Simulator is a kind of Separated SpO2 simulator, small and light. It can perform a series of tests for the oximeter by simulation means, and gives cognizance of veracity about the oximeter. Because the oximeter manufacturer may use different R-curves, some R-curves used abroad currently have been embedded.

Features

- TFT-LCD With Adjustable Backlight Brightness.
- Separated design for simulator probe and host.
- Film Button.
- Simulate presettable patient state.
- Test Reaction Time of The Device.
- Simulate SPO2 and PR Under Different Amplitudes.
- Test Device Performance Under Different Interference Source.
- With Rechargeable Lithium Battery.