

NIBP Simulator



Model : YR02282
NIBP Simulator
Manufacturer : Kalstein

Price : \$ 6,000.00

Description :

Introduction

NIBP Simulator is a multi-purpose test instrument for use with oscillometric Non-Invasive Blood Pressure(NIBP) Monitors. The appearance is concise, operation is simple and convenient, with its internal pump, the Simulator can generate pressures up to 400 mmHg (53.3 kPa). The Simulator provides dynamic blood pressure simulations, static pressure calibration, automated leak testing, and pressure relief valve testing. You can quickly and exactly analyse different NIBP Monitors.

Features

- 4.3"LCD, Test Functions for Pressure Gauge, Pressure Leak and Relief Valve.

- Built-in Pump Provides Pressure Source.

- Simulate Ambulatory Blood Pressure for Arm and Wrist Sphygmomanometer.

- Simulator 8 Kinds of Patient State and User-defined.

- Built-in Air Chamber to Simulate Adult and Neonate Cuff.