

Constant Temperature Incubator



Model : YR02044
Constant Temperature Incubator
Manufacturer : Kalstein
Price : \$ 0.00
Description :

Thermostatic Incubator

Large integrated observation window; fan circulation blast, reduce temperature deviation and fluctuation; adopt magnetic door seal, keep sealing in the box; shelf spacing, quantity adjustable.

Suitable for the use of routine constant temperature culture.

Use for constant temperature culture in the medical and health, pharmaceutical industry, biochemistry and agricultural research.

Product introduction?

For medical health, pharmaceutical, chemical and agricultural and other departments as constant temperature incubator.

Feature?

-
Easy to watch the goods in the inner

chamber

with integrated window.

-

Constant Temperature Incubator

Equipped with circulation fan in the inner chamber to reduce the temperature deviation and volatility.

- Magnetic door seals to keep the liner chamber sealed.

- Shelf can be adjusted.

Parameters by Models

Constant Temperature Incubator

Model	YR02044	YR02045	YR02046	YR02047
Temp Range (° C)	RT+5~65			
Accuracy (° C)	±1			
Fluctuation (%)	± 1			
Voltage (V)	220			
Power (KW)	0.3			
Inner Size(mm)	250×250×250	420×360×360	500×400×400	600×500×500
Overall Dimension (mm)	480×520×400	610×580×470	690×620×510	790×720×610
Packing Size (mm)	600×600×480	730×660×550	810×700×590	910×800×690
Net Weight (KG)	16.5	23.5	29.5	41.5
Shipping Weight (KG)	18.5	26.5	33.5	44.5

Constant Temperature Incubator

Options available for Constant Temperature Incubator :

[check price](#)

(YR02045) 54 Liters (+\$ 93.00),
(YR02046) 80 Liters (+\$ 155.00),
(YR02047) 150 Liters (+\$ 341.00),
(YR02044) 15 Liters.