

Trinocular Microscope



Model : YR0245
Trinocular Microscope
Manufacturer : -
Price : \$ 5,920.00
Description :

Trinocular Microscope

Description:

Optical System Infinite Optical System.

Viewing Head Seidentopf Binocular Viewing Head, Inclined at 45°, 360° Rotatable, Interpupillary 48-75mm.

Eyepiece

?30 Wide Field Eyepiece WF10×/22, Eyepiece Tube Diameter 30mm.

Wide Field Eyepiece WF15×/16

Wide Field Eyepiece WF20×/12

Parameters by Model

Trinocular Microscope

Model	YR0245		
Optical System	Infinite Optical System	●	
Viewing Head	Top Binocular Viewing Head, Inclined at 45°, 360° Rotatable, Interpupillary 48	●	
Eyepiece	Φ30 Wide Field Eyepiece WF10×/22, Eyepiece Tube Diameter 30mm	●	
	Wide Field Eyepiece WF15×/16	○	
	Wide Field Eyepiece WF20×/12	○	
Objective	Infinite Plan Achromatic Objective 4×/0.10 WD=22	●	
	LWD Infinite Plan Achromatic Phase Objective	10×/0.25 WD=6	●
		20×/0.4 WD=3.1	●
		40×/0.55 WD=2.2	●
	Lamp House Adjustment Objective	○	
Nosepiece	Backward Quintuple Nosepiece	●	
Condenser	ELWD Condenser NA0.3, LWD72mm (without Condenser 150mm)	●	
Centering Telescope	Centering Telescope (φ30mm)	●	
Phase Annulus	10×, 20×, 40× Phase Annulus Plate (Adjustable)	●	
Stage	Plain Stage 230×170	●	
	Glasse Insert	●	
	Attachable Mechanical Stage, X, Y Coaxial Control, Moving Rang 80×120mm	●	
	Auxiliary Stage	●	
	Terasaki Holder	●	
	Petri Dish Holder φ38	●	
	Petri Dish Holder φ54	●	
Focusing System	Coaxial Coarse and Fine Adjustment, Fine Division 0.002mm, Moving Range up 4.5mm down 4.5mm	●	
Illumination	6V30W Halogen Lamp, Brightness Adjustable	●	
Filter	Blue, Green and Ground Glasse, Diameter 45mm	●	
Accessories	Photo Attachment	○	
	C Mount 0.4×	○	
	Epi-Fluorescent Attachment	○	

Note: ● Standard Outfit, ○ Optional