

## ROM Pneumatic Walker



**Model : YR02411**

**ROM Pneumatic Walker**

**Manufacturer : Kalstein**

**Price : \$ 144.00**

**Description :**

ROM Pneumatic Walker

Indications:

- Grade II/? ankle sprain
- Stable fracture of the fibula, talus, calcaneus and medial malleolus
- Post ligament, soft tissue and tendon surgery
- Forefoot or midfoot injury

Features:

- Rocker bottom design assists gait
- Dial lock hinges with adjustable ROM limitation
- Dorsiflexion limitation at: 0°, 7.5°, 15°, 22.5°, 30°, 37.5°, 45°
- Plantarflexion limitation at: 0°, 7.5°, 15°, 22.5°, 30°, 37.5°, 45°
- Immobilization limitation at: 0°, 7.5°, 15°, 22.5°, 30°, 37.5°, 45°
- Available in 14" & 17" height
- Inflatable aircells are available for further support while reducing swelling.

Size: XS, S, M, L, XL