

ROM Elbow Brace with Sling



Model : YR02479
ROM Elbow Brace with Sling
Manufacturer : Kalstein

Price : \$ 97.50

Description :

ROM Elbow Brace with Sling

Indications:

- ? Post-operative immobilization or ROM control for the elbow
- ? Conservative treatment of elbow dislocation or luxation
- ? Stable or internally fixed fractures of the distal humerus or proximal radius or ulna

Features:

- ? Immobilization or fix the elbow in every 15°
- ? With padded shoulder strap to release the pressure around shoulder
- ? Simply press and rotate hinge with lock
- ? Extension limitation from 0° ~ 60°
- ? Flexion limitation from 0° ~ 120°
- ? Flexion from 0° ~ 60°
- ? Adjustable in every 15°

Size: Universal (Left/Right)